

BONE SHARK[®]

Disposable Bone Mill

Engineered for efficiency, safety, and power — with zero cleanup.
Single-use. Powerful. Ready to perform.

Reimagining Bone Milling

The Bone Shark[®] Disposable Bone Mill is the healthcare industry's leading single-use bone milling device — designed to deliver sterile, safe, and powerful performance with no need for reprocessing equipment.

NEW Patented design highlight:

- Larger loading chamber and opening making it faster & easier to load & process bone material
- Double the diameter, re-engineered end mill cutter with ability to cut in both directions. Delivers faster, smoother, and more efficient milling with every pass
- New optional integration with standard OR drill — you can now power the system, streamlining your workflow.



Scan the QR Code to
View the Product Video


LEEmedical
Advancing Innovation



The Smarter Way to Mill Autograft Bone
Single-use. Powerful. Ready to perform.

BONESHARK® Disposable Bone Mill



Scan the QR Code to
View the Product Video

LEEmedical
Advancing Innovation

Bone Shark is a Registered Trademark of Tarca LLC.
Patent Numbers: US7156329, 11701243
© Copyright 2025 Lee Medical. All Rights Reserved.



Smarter Bone Harvesting, Single-Use Simplicity

Key Benefits at a Glance

- ☒ 57% Larger Chamber Capacity – Accepts fragments up to 1.5 cm for increased throughput.
- ☒ Dual-Direction Cutting – A 100% larger 4-fluted mill cuts in both directions for efficiency.
- ☒ Drill Integration – Use manually or connect to a standard OR drill in oscillate mode.
- ☒ Pre-Assembled & Sterile – Open, use, and dispose with no reprocessing required.

Ideal Clinical Applications

- ☒ Spinal Fusion
- ☒ Orthopedic Reconstruction
- ☒ ENT Procedures
- ☒ Dental Implants
- ☒ Maxillofacial Surgery
- ☒ Small Bone Applications

For more information please view our video on www.leemedicalUS.com or call Customer Service at 888-422-7717

 888.422.7717

 www.leemedicalUS.com

 info@leemedicalnj.com

 30 Vreeland Drive, Suite 8
Skillman, NJ, 08558